

# Bacterial Genomic DNA Isolation Kit

rev 07/20

(Catalog # K309-100; 100 isolations; Store at Multiple Temperatures)

- I. Introduction:**  
 Bacteria are one of the most abundant and diverse organisms on the planet, which take part in numerous critical ecosystem processes. Many bacterial species are pathogens that are responsible for causing a variety of human and animal diseases. In addition to their medical and ecological importance, bacteria are also used in various industrial applications such as production of enzymes and biofuels. BioVision's bacterial genomic DNA isolation kit provides convenient and simple step-by-step method for isolating quality genomic DNA from gram-negative and gram-positive bacterial species. This kit utilizes enzymatic reactions to release bacterial DNA from the cell. DNA release from the cell is coupled with adsorption of DNA onto a silica spin-column in the presence of high salt concentration, eliminating the use of toxic organic compounds or solvents. DNA purified by this kit is suitable for various downstream molecular biology applications such as PCR, cloning, DNA hybridization, and Southern Blotting.
- II. Applications:**
- PCR
  - Cloning
  - DNA Hybridization
  - Southern Blotting
- III. Sample Type:**
- Gram positive and gram negative bacterial species
- IV. Kit Contents:**

| Components                      | K309-100  | Cap Code | Part Number | Storage (°C) |
|---------------------------------|-----------|----------|-------------|--------------|
| Buffer A [Re-suspension Buffer] | 25 ml     | NM/Clear | K309-100-1  | RT           |
| Enzyme Mix A                    | 1 ml      | Green    | K309-100-2  | -20°C        |
| RNase A                         | 600 µl    | Blue     | K309-100-3  | -20°C        |
| Buffer B [Reaction Buffer]      | 1.7 ml    | Yellow   | K309-100-4  | RT           |
| Enzyme Mix B                    | 1.2 ml    | Red      | K309-100-5  | -20°C        |
| Buffer C [Binding Buffer]*      | 25 ml     | NM/Brown | K309-100-6  | RT           |
| Buffer D [Wash Buffer]**        | 30 ml     | WM       | K309-100-7  | RT           |
| Buffer E [Elution Buffer]       | 22 ml     | WM       | K309-100-8  | RT           |
| Spin Columns/Collection Tubes   | 100 tubes | -        | K309-100-9  | RT           |

\*Add 28 mL 100% Ethanol to the Buffer C [Binding Buffer] bottle and

\*\*136 mL 100% Ethanol to the Buffer D [Wash Buffer] bottle before use.

- V. User Supplied Reagents and Equipment:**
- DNase-free aerosol tips and micro-centrifuge tubes
  - 100% Ethanol
  - Heating Block
  - Centrifuge
  - Mutanolysin and Lysostaphin
- VI. Storage Conditions and Reagent Preparation:**  
 Store kit at -20°C and room temperature (RT), protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.
- **Buffer A, B, and E:** Ready to use. Store at RT.
  - **Buffer C:** Add 28 ml of 100% Ethanol, molecular biology grade. Mix well and store at RT.
  - **Buffer D:** Add 136 ml of 100% Ethanol, molecular biology grade. Mix well and store at RT.
  - **Enzyme Mix A, Enzyme Mix B, and RNaseA:** Ready to use. Store at -20°C. Keep on ice at all times while in use.
  - **Spin Columns:** Ready to use. Store at RT in dry conditions.
- VII. Bacterial Genomic DNA Extraction Protocol:**
- 1. Sample Preparation:**
- Transfer 1ml of overnight bacterial culture into a 1.5 ml tube and pellet the cells by centrifugation at 12,000 x g for 1 min at 4°C. Discard supernatant.
- Note:** If using >1-1.5 ml overnight culture, scale up the entire protocol proportionally.
- Re-suspend the pellet in 250 µl Buffer A [Re-suspension Buffer].




Figure 1.

1% 1X TBE agarose gel

**Figure 1. Example genomic DNA isolation using BioVision's Bacterial Genomic DNA Isolation Kit.** BioVision's Bacterial Genomic DNA Isolation kit efficiently purifies high quality genomic DNA from gram negative bacteria (such as *Escherichia coli*) and gram positive bacteria (such as *Staphylococcus aureus* and *Micrococcus luteus*). A: BriteRuler 1kb DNA Ladder; B: *M. luteus*; C: *S. aureus*; D: *E. coli*.

- c. Add 10  $\mu$ l of Enzyme Mix A. If RNA free DNA is desired, add 6  $\mu$ l of RNase A solution at this time. Mix by inverting the tube 3-5 times.

**Note:** Certain bacterial species require additional enzymes during cell lysis. For *Streptococcus mutans* and *Staphylococcus species*, supplement the lysis reaction with  $\geq 250$  units/mL Mutanolysin\* and  $\geq 200$  units/mL Lysostaphin\* respectively. \*Not supplied.

- d. Incubate for 30 min at 37°C.

## 2. DNA Release:

- e. Add 17  $\mu$ l of Buffer B [Reaction Buffer] and 12  $\mu$ l of Enzyme Mix B. Mix by inverting the tube 3-5 times.  
 f. Incubate for 30 min at 55°C.

## 3. Binding:

- g. Add 500  $\mu$ l of Buffer C [Binding Buffer] and mix by pipetting up and down until the solution becomes clear.  
 h. Place the provided spin column into the provided collection tube and pipette the entire supernatant onto the top of the column.  
 i. Spin at 12,000 x g for 1 min at 4°C and discard the flow through.

## 4. Washing:

- j. Add 750  $\mu$ l of Buffer D [Wash Buffer] onto the top of the spin column and spin (12,000 x g; 1 min., 4°C). Discard the flow through.  
 k. Repeat step j one more time.  
 l. Spin the spin column at 12,000 x g for 2 min at 4°C to dry.

## 5. Elution:


- m. Transfer the spin column to a clean, DNase-free 1.5ml tube.  
 n. Add 200  $\mu$ l of Buffer E [Elution Buffer] to the top of the spin column and incubate for 1-2 min at RT.  
 o. Spin at 12,000 x g for 1 min at 4°C to elute the genomic DNA.  
 p. Sample is now purified and ready to use. Store genomic DNA at -20°C or immediately use the sample in a downstream application of your choice.

**Note:** \*(Note: Generally, good quality genomic DNA will have A260/280 of 1.7- 1.99 and exhibit one clear band of high molecular weight on 1% agarose gel. See troubleshooting guide in Section VIII for help

## VIII. Troubleshooting:

| Issue                           | Possible Reason            | Recommendations  |
|---------------------------------|----------------------------|--|
| Low yield                       | Low bacteria concentration | Monitor OD <sub>600 nm</sub> of your overnight culture. We recommend OD <sub>600 nm</sub> range to be between 1 and 2. |
|                                 | Incomplete Lysis           | Increase incubation with Enzyme Mix A up to 45 min - 1hr.  |
|                                 | Incomplete DNA release     | Increase incubation with Enzyme Mix B up to 45 min - 1hr.  |
| Low A260/280 (<1.6)             | Protein Contamination      | Increase incubation with Enzyme Mix B up to 45 min - 1hr.  |
| High A260/280 (>2.0)            | RNA Contamination          | Add RNase A during the cell lysis step.  |
| No DNA band/severe smear on gel | DNase contamination        | Use DNase free aerosol tips, DNase-free tubes, and practice good sterile technique                                     |

## Protocol Quick Guide



**FOR RESEARCH USE ONLY!**  
**Not to be used on humans**